

DRUG PARENT CHART WITH CUT-OFF LEVELS

PANEL STRIPS Blue bars below are test strip panels built into Identify® Diagnostics cup or dip tests. Your test may include one or more panels.

BRAND® / STREET NAMES Green tabs below are common drug brands or street names. Not all drug brands or street names are included in this list.

Drug Compound	Cut-off in ng/mL	Drug Compound	Cut-off in ng/mL
AMPHETAMINE 1000 (AMP)		METHADONE 300 (MTD)	
D-Amphetamine ADDERALL®	1000	Methadone DOLOPHINE® METHADOSE®	300
D,LAmphetamine	2500	(±)2-Ethyl-1,5-dimethyl-3,3-diphenylpyrrolinium	50000
L-Amphetamine	50000	Doxylamine DICLECTIN®	50000
(±)3,4-Methylenedioxyamphetamine (MDA)	2000	METHAMPHETAMINES 1000 (mAMP/MET)	
Ephedrine	>100000	+/-Methamphetamine DESOXYN®	2000
3,4-Methylenedioxythylamphetamine (MDEA)	>100000	+Methamphetamine	1000
BARBITURATES 300 (BAR)		3,4-Methylenedioxyethylamphetamine (MDEA)	35000
Secobarbital SECONAL®	300	(+/-)3,4-Methylenedioxymethamphetamine (MDMA)	2000
Phenobarbital LUMINAL®	2500	Ranitidine (Zantac)	>100000
Butalbital FORTABS®	500	3,4-Methylenedioxyamphetamine (MDA)	>100000
Pentobarbital NEMBUTAL®	1500	D-Amphetamine	>100000
Amobarbital AMYTAL®	2500	L-Amphetamine	>100000
Cyclopentobarbital	500	Ephedrine	>100000
Butethal FIORICET®	800	METHYLENEDIOXYMETHAMPHETAMINE 500 (MDMA)	
Barbital VEVONAL®	300	(+/-)3,4-Methylenedioxymethamphetamine (MDMA) ECSTASY	500
Butabarbital BUTISOL®	1500	D-Amphetamine	>100000
BENZODIAZEPINES 300 (BZO)		L-Methamphetamine	100000
Oxazepam SERAX®	300	3,4-Methylenedioxyethylamphetamine (MDEA)	200
Alprazolam XANAX®	200	3,4-Methylenedioxyamphetamine (MDA)	2000
α-Hydroxyalprazolam	1000	OPIATE 300 (MOP)	
Bromazepam LEXOTANIL® LECTOPAM®	250	Morphine HEROIN METABOLITE ROXANOL®	300
Chlordiazepoxide LIBRIUM®	2500	Codeine TYLENOL 3®	300
Clobazam ONFI®	100	Hydrocodone LORATAB® LORCET® VICODIN®	2000
Clonazepam KLONOPIN®	850	Hydromorphone DILAUDID® OPIDOL®	3500
Clorazepate TRANXENE®	250	Morphine 3-β-D-glucuronide	300
Delorazepam DADUMIR®	1600	6-Monoacetylmorphine HEROIN METABOLITE	600
Diazepam VALIUM®	200	Normorphone	100000
Estazolam PROSOM®	200	Oxycodone OXYCONTIN® PERCOCET®	10000
Flunitrazepam	300	Oxymorphone OPANA®	50000
Lorazepam ATIVAN®	1000	Thebaine	7000
Midazolam	1500	OXYCODONE 100 (OXY)	
Nitrazepam	100	Oxycodone OXYCONTIN® PERCOCET®	100
Nordiazepam	400	Morphine HEROIN METABOLITE ROXANOL®	50000
Temazepam RESTORIL®	150	Codeine TYLENOL 3®	25000
Triazolam HALCION®	500	Morphine 3-β-D-glucuronide	50000
BUPRENORPHINE 10 (BUP)		Hydrocodone LORATAB® LORCET® VICODIN®	1600
Buprenorphine SUBOXONE® SUBUTEX®	10	Hydromorphone DILAUDID® OPIDOL®	15000
Norbuprenorphine	15	Normorphone	100000
Buprenorphine-3-D-glucuronide	12.5	Oxymorphone OPANA®	1500
Norbuprenorphine-3-D-glucuronide	175	PHENCYCLIDINE 25 (PCP)	
Morphine-3-D-glucuronide	100000	Phencyclidine ANGEL DUST	25
Morphine HEROIN METABOLITE	>100000	4-Hydroxyphencyclidine	15000
Oxymorphone OPANA®	>100000	TRICYCLIC ANTIDEPRESSANTS 1000 (TCA)	
Hydromorphone DILAUDID® OPIDOL®	>100000	Notriptiline PAMELOR®	1000
COCAINE 300 (COC)		Trimipramine SURMONTIL®	4500
Cocaine COKE, CRACK	>100000	Amitriptyline ELAVIL®	1000
Benzoylecgonine	300	Promazine SPARINE®	3000
Ecgonine HCl	35000	Desipramine NORPRAMIN®	1000
MARIJUANA 50 (THC)		Imipramine TOFRANIL®	1000
11-nor-Δ ⁸ -THC-9-COOH POT, WEED	50	Clomipramine ANAFRANIL®	7500
11-nor-Δ ⁹ -THC-9-COOH POT, WEED	50	Doxepin SILENOR® ZONALON®	3000
Δ ⁸ -Tetrahydrocannabinol	8000	Maprotiline LUDIOMIL®	50000
Δ ⁹ -Tetrahydrocannabinol	10000		
Cannabinol	10000		
Cannabidiol	100000		

The data presented in this chart demonstrates excellent precision of Identify® Diagnostics tests across multiple concentrations of human urine. Analytical Specificity: Cross-reactivity was established by spiking various concentrations of similarly structured drug compounds into drug-free urine / a negative control. Analyzing concentrations of each compound by using Identify® Diagnostics tests, the concentration of the drug that produced a response approximately equivalent to the cut-off concentration of the assay was determined. This is not a CLIA document.