



# StrepAim®

A rapid DipStick test for the Qualitative Detection of Group A Streptococcal Antigen from throat swabs.

COMPARES WITH NATIONALLY RECOGNISED BRANDS.

**Simple, Reliable  
Strep A Results!**

LOOK AT THESE GREAT FEATURES!



**CLIA  
WAIVED**



**99% ACCURACY**

Provides accurate results that are easy to read.



Convenient workstation, dipsticks, disposable test tubes, sterile swabs, positive and negative controls.



Results in 5 minutes.



Store at 2° to 30°C.

**StrepAim® KIT SIZE**



	SIZE	REF#
<b>StrepAim®</b>	25 Tests	73025



➔ CAP / CLSI formatted procedures available via fax or email - just add your name!

CPT CODE: 87880 QW..... Reimbursement \$16.76\*



**Medical Distribution Group, Inc.**

For questions please call: 727-744-2967

Medical Distribution Group, Inc. 2361 Whitfield Park Avenue. Sarasota, FL 34243