

DRUG PARENT CHART

PANEL STRIPS

Blue bars below represent test panels strips built into Identify Diagnostics cup and dip tests. Your test may include one or more panels based on your specific configuration.

BRAND® / STREET NAMES

Green tabs below are common drug brands or street names. Not all drug brands or street names are included in this list. The list of drug compounds under each blue bar indicate what drugs are detectable by that panel and at what cut-off level in ng/mL.

DRUG COMPOUND	Cut-off in ng/ml	DRUG COMPOUND	Cut-off in ng/ml
AMPHETAMINE 1000 AMP		MARIJUANA 50 THC	
D-Amphetamine ADDERALL®	1,000	11-nor- Δ^8 -THC-9-COOH POT, WEED	50
D,LAmphetamine	2,500	11-nor- Δ^9 -THC-9-COOH POT, WEED	50
L-Amphetamine	50,000	Δ^8 -Tetrahydrocannabinol	8,000
(\pm)3,4-Methylenedioxyamphetamine (MDA)	2,000	Δ^9 -Tetrahydrocannabinol	10,000
Ephedrine	>100,000	Cannabinol	10,000
3,4-Methylenedioxyethylamphetamine (MDEA)	>100,000	Cannabidiol	100,000
BARBITURATES 300 BAR		METHADONE 300 MTD	
Secobarbital SECONAL®	300	Methadone DOLOPHINE® METHADOSE®	300
Phenobarbital LUMINAL®	2,500	Doxylamine DICLECTIN®	50,000
Butalbital FORTABS®	500	METHAMPHETAMINES 1000 MET	
Pentobarbital NEMBUTAL®	1,500	+/-Methamphetamine DESOXYN®	2,000
Amobarbital AMYTAL®	2,500	+Methamphetamine	1,000
Cyclopentobarbital	500	3,4-Methylenedioxyethylamphetamine	35,000
Butethal FIORICET®	800	(+/-)3,4-Methylenedioxymethamphetamine (MDMA)	2,000
Barbital VEVONAL®	300	Ranitidine ZANTAC®	>100,000
Butobarbital BUTISOL®	1,500	3,4-Methylenedioxyamphetamine (MDA)	>100,000
BENZODIAZEPINES 300 BZO		D-Amphetamine	>100,000
Oxazepam SERAX®	300	L-Amphetamine	>100,000
Alprazolam XANAX®	200	Ephedrine	>100,000
a-Hydroxyalprazolam	1,000	METHYLENEDIOXYMETHAMPHETAMINE 500 MDMA	
Bromazepam LEXOTANIL® LECTOPAM®	250	(+/-)3,4-Methylenedioxymethamphetamine (MDMA) ECSTASY	500
Chlordiazepoxide LIBRIUM®	2,500	D-Amphetamine	>100,000
Clobazam ONFI®	100	L-Methamphetamine	100,000
Clonazepam KLONOPIN®	850	3,4-Methylenedioxyethylamphetamine (MDEA)	200
Clorazepate TRANXENE®	250	3,4-Methylenedioxyamphetamine (MDA)	2,000
Delorazepam DADUMIR®	1,600	OPIATE 300 MOP	
Diazepam VALIUM®	200	Morphine HEROIN METABOLITE ROXANOL®	300
Estazolam PROSOM®	200	Codeine TYLENOL 3®	300
Flunitrazepam	300	Hydrocodone LORATAB® LORCET® VICODIN®	2,000
Lorazepam ATIVAN®	1,000	Hydromorphone DILAUDID®	3,500
Midazolam	1,500	Morphine 3- β -D-glucuronide	300
Nitrazepam	100	6-Monoacetylmorphine HEROIN METABOLITE	600
Nordiazepam	400	Normorphine	100,000
Temazepam RESTORIL®	150	Oxycodone OXYCONTIN®	10,000
Triazolam HALCION®	500	Oxymorphone OPANA®	50,000
BUPRENORPHINE 10 BUP		Thebaine	7,000
Buprenorphine SUBOXONE® SUBUTEX®	10	OXYCODONE 100 OXY	
Norbuprenorphine	15	Oxycodone OXYCONTIN® PERCOCET®	100
Buprenorphine-3-D-glucuronide	12.5	Morpine HEROIN METABOLITE ROXANOL®	50,000
Norpuprenorphine-3-D-glucuronide	175	Codeine TYLENOL 3®	25,000
Morphine-3-D-glucuronide	100,000	Morphine 3- β -D-glucuronide	50,000
Morpine HEROIN METABOLITE	>100,000	Hydrocodone LORATAB® LORCET® VICODIN®	1,600
Oxymorphone OPANA®	>100,000	Hydromorphone DILAUDID®	15,000
Hydromorphone DILAUDID® OPIDOL®	>100,000	Normorphine	100,000
COCAINE 300 COC		Oxymorphone OPANA®	1,500
Cocaine COKE, CRACK	>100,000	PHENCYCLIDINE 25 PCP	
Benzoylcegonine COCAINE METABOLITE	300	Phencyclidine ANGEL DUST	25
Ecgonine HCl	35,000	SYNTHETIC MARIJUANA 25 K2	
COTININE 200 COT		Synthetic Marijuana SPICE	25
Cotinine NICOTINE	200	TRAMADOL 100 TRA	
EDDP 300 EDDP		Tramadol CONZIP® ULTRACET® ULTRAM®	100
EDDP Methadone Metabolite	300	TRICYCLIC ANTIDEPRESSANTS 1000 TCA	
ETHYL GLUCURONIDE 500 ETG		Notriptiline PAMELOR®	1,000
Ethyl Glucuronide ALCOHOL	500	Trimipramine SURMONTIL®	4,500
FENTANYL 25 FEN		Amitriptyline ELAVIL®	1,000
Fentanyl DURAGESIC® ABSTRAL® LONSY®	25	Promazine SPARINE®	3,000
GABAPENTIN 2000 GAB		Desipramine NORPRAMIN®	1,000
Gabapentin GRALISE® NEURONTIN®	2,000	Imipramine TOFRANIL®	1,000
KRATOM 100 KRA		Clomipramine ANAFRANIL®	7,500
Kratom	100	Doxepin SILENOR® ZONALON®	3,000
		Maprotiline LUDIOMIL®	50,000