

# DRUG PARENT CHART

## PANEL STRIPS

Blue bars below represent test panels strips built into Identify **HEALTH** cup and dip tests. Your test may include one or more panels based on your specific configuration.

## BRAND® / STREET NAMES

Green tabs below are common drug brands or street names. Not all drug brands or street names are included in this list. The list of drug compounds under each blue bar indicate what drugs are detectable by that panel and at what cut-off level in ng/mL.

DRUG COMPOUND	Cut-off in ng/ml	DRUG COMPOUND	Cut-off in ng/ml
<b>AMPHETAMINE 500 AMP</b>		<b>MARIJUANA 50 THC</b>	
D-Amphetamine <b>ADDERALL®</b>	500	11-nor- $\Delta^9$ -THC-9-COOH <b>POT, WEED</b>	50
L-Amphetamine	50,000	(+/-)-11-Hydroxy- $\Delta^9$ -THC	5,000
Phentermine	45,000	(-)- $\Delta^8$ -THC	20,000
MDA	8,000	(-)- $\Delta^9$ -THC	20,000
<b>BARBITURATES 300 BAR</b>		<b>METHADONE 300 MTD</b>	
Secobarbital <b>SECONAL®</b>	300	Methadone <b>DOLOPHINE®</b> <b>METHADOSE®</b>	300
Phenobarbital <b>LUMINAL®</b>	300		
Butalbital <b>FORTABS®</b>	300	<b>METHAMPHETAMINES 500 MET</b>	
Pentobarbital <b>NEMBUTAL®</b>	250	D-Methamphetamine <b>DESOXYN®</b>	500
Amobarbital <b>AMYTAL®</b>	2,500	D-Methamphetamine	50,000
Cyclopentobarbital	500	L-Amphetamine	50,000
Aprobarbital	500	1R,2S(-)-Ephedrine	100,000
Butabarbital <b>BUTISOL®</b>	100	MDMA	3,500
<b>BENZODIAZEPINES 300 BZO</b>		<b>METHYLENEDIOXYMETHAMPHETAMINE 500 MDMA</b>	
Oxazepam <b>SERAX®</b>	300	(+/-)-MDMA <b>ECSTASY</b>	500
Alprazolam <b>XANAX®</b>	200	(+/-)-MDA	3,900
a-Hydroxyalprazolam	1,900	(+/-)-MDEA	500
Bromazepam <b>LECTOPAM®</b> <b>LEXOTANIL®</b>	1,000		
Clobazam <b>ONFI®</b>	200	<b>OPIATE 300 OPI300</b>	
Clorazepate <b>TRANXENE®</b>	750	Morphine <b>HEROIN METABOLITE</b> <b>ROXANOL®</b>	300
Diazepam <b>VALIUM®</b>	1,000	Codeine <b>TYLENOL 3®</b>	100
Flunitrazepam	250	Hydrocodone <b>LORATAB®</b> <b>LORCET®</b> <b>VICODIN®</b>	1,250
Lorazepam <b>ATIVAN®</b>	3,900	Hydromorphone <b>DILAUDID®</b> <b>OPIDOL®</b>	2,500
Lorazepam-Glucuronide	5,000	Levorphanol	50,000
Nitrazepam	250	Morphine 3-Glucuronide	400
Norchlordiazepoxide	500	Norcodeine	6,000
Nordazepam	390	Oxycodone <b>OXYCONTIN®</b> <b>PERCOCET®</b>	75,000
Nordiazepoxide	400	Thebaine	90,000
Temazepam <b>RESTORIL®</b>	150		
Triazolam <b>HALCION®</b>	2,500	<b>OXYCODONE 100 OXY</b>	
		Oxycodone <b>OXYCONTIN®</b> <b>PERCOCET®</b>	100
<b>BUPRENORPHINE 10 BUP</b>		Codeine <b>TYLENOL 3®</b>	50,000
Buprenorphine <b>SUBOXONE®</b> <b>SUBUTEX®</b>	10	Ethylmorphine	50,000
		Hydrocodone <b>LORATAB®</b> <b>LORCET®</b> <b>VICODIN®</b>	5,000
<b>COCAINE 150 COC</b>		Hydromorphone <b>DILAUDID®</b> <b>OPIDOL®</b>	25,000
Cocaine <b>COKE, CRACK</b>	5,000	Oxymorphone <b>OPANA®</b>	12,500
Benzoylcegonine	150		
Ecgonine	50,000	<b>PHENCYCLIDINE 25 PCP</b>	
Cocaethylene	50,000	Phencyclidine <b>ANGEL DUST</b>	25
<b>COTININE 200 COT</b>		<b>TRAMADOL 100 TRA</b>	
Cotinine <b>NICOTINE</b>	200	Tramadol <b>CONZIP®</b> <b>ULTRACET®</b> <b>ULTRAM®</b>	100
<b>EDDP 300 EDDP</b>		<b>SYNTHETIC MARIJUANA 25 K2</b>	
EDDP Methadone Metabolite	300	Synthetic Marijuana <b>SPICE</b>	25
<b>ETHYL CLUCURONIDE 500 ETG</b>		<b>TRICYCLIC ANTIDEPRESSANTS 1000 TCA</b>	
Ethyl Glucuronide <b>ALCOHOL</b>	500	Notriptyline <b>PAMELOR®</b>	1,000
		Amitriptyline <b>ELAVIL®</b>	4,000
<b>FENTANYL 25 FEN</b>		Clomipramine <b>ANAFRANIL®</b>	2,000
Fentanyl <b>DURAGESIC®</b> <b>ABSTRAL®</b> <b>LONSY®</b>	25	Desipramine <b>NORPRAMIN®</b>	500
		Doxepin <b>SILENOR®</b> <b>ZONALON®</b>	1,000
<b>GABAPENTIN 2000 GAB</b>		Imipramine <b>TOFRANIL®</b>	1,000
Gabapentin <b>GRALISE®</b> <b>NEURONTIN®</b>	2,000	Promethazine	1,000
		Trimipramine <b>SURMONTIL®</b>	5,000
<b>KRATOM 100 KRA</b>			
Kratom	100		