

ALL CUP TESTS

WITH CUT-OFF LEVELS

ID-CP6  ID-US6  ID-US6-BUP  ID-H7-1 ID-H9-1 ID-CP10 ID-CP12-BUP ID-CP12-BUP-ADULT  ID-US12-BUP  ID-H12-EFKK-1 ID-H12-1 ID-H12-2 ID-H13-1 ID-US14-ADULT  ID-H14-ADULT-EFKT-2 ID-H16-ADULT-EFKT ID-US16-ADULT-EFKT 


SUBSTANCE TESTED	CUP TESTS																
	ID-CP6	ID-US6	ID-US6-BUP	ID-H7-1	ID-H9-1	ID-CP10	ID-CP12-BUP	ID-CP12-BUP-ADULT	ID-US12-BUP	ID-H12-EFKK-1	ID-H12-1	ID-H12-2	ID-H13-1	ID-US14-ADULT	ID-H14-ADULT-EFKT-2	ID-H16-ADULT-EFKT	ID-US16-ADULT-EFKT
Amphetamines (AMP)	1000	500		500	500	1000	1000	1000	500	500	500	500	500	500	500	500	500
Barbiturates (BAR)					300	300	300	300	300		300	300	300	300		300	300
Benzodiazepines (BZO)	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300
Buprenorphine (BUP)			10	10	10		10	10	10	10	10	10	10	10	10	10	10
Cocaine (COC)	300	150	150	150	150	300	300	300	150	150	150	150	150	150	150	150	150
Ecstasy (MDMA)						500	500	500		500	500	500	500	500	500	500	500
Marijuana (THC)	50	50		50		50	50	50	50	50	50	50	50	50	50	50	50
Methadone (MTD)			300		300	300	300	300	300		300	300	300	300	300		
Methadone Metabolite (EDDP)														300		300	300
Methamphetamines (MET)					500	1000	1000	1000	500	500	500	500	500	500	500	500	500
Opiates (MOP) or (OPI300)	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300
Oxycodone (OXY)	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100
Phencyclidine (PCP)						25	25	25	25		25		25	25		25	25
Tricyclic Antidepressants (TCA)												1000	1000	1000			
** Cotinine / Nicotine (COT)																	
** Ethyl Glucuronide (ETG)										500					500	500	500
** Fentanyl (FEN)										20					20	25	25
** Gabapentin (GAB)																	
** Kratom (KRA)										100							
** Synthetic Marijuana (K2)										25					25	25	25
** Tramadol (TRA)															100	100	100
*ADULTERANTS																	
Creatinine (C)								✓						✓	✓	✓	✓
Glutaraldehyde (G)														✓		✓	✓
Acidic or Alkaline pH (P)								✓						✓	✓	✓	✓
Nitrite (N)														✓		✓	✓
Oxidants (O)														✓	✓	✓	✓
Specific Gravity (S)								✓						✓		✓	✓

** NOTE: Premium panels not CLIA Waived intended for forensic use only. * NOTE: Urine tests with adulterants measure levels of substances in human urine which may indicate sample tampering.

A green box  indicates the substance or adulterant in the far left column is tested for. The number inside the box indicates the cut-off level in ng/ml. A cut-off level is the concentration a med or narcotic must exceed in the urine sample to trigger a positive result.

Tests above with the American Flag  are designed & manufactured in the United States. All Identify Diagnostics & Identify Diagnostics USA tests are made with USA reagents. Reagents are the substances on each panel strip which produce the urine test result.

ALL DIP TESTS WITH CUT-OFF LEVELS



SUBSTANCE TESTED	ID-US1-DIP-COT	ID-US1-DIP-THC	ID-US1-DIP-ETG	ID-US1-DIP-FEN	ID-US1-DIP-GAB	ID-US1-DIP-KRA	ID-CP5-DIP	ID-CP6-DIP	ID-CP10-DIP	ID-US10-DIP-1	ID-CP12-DIP
Amphetamines (AMP)						1000		1000	500	1000	
Barbiturates (BAR)								300	300	300	
Benzodiazepines (BZO)							300	300	300	300	
Buprenorphine (BUP)									10	10	
Cocaine (COC)						300	300	300	150	300	
Ecstasy (MDMA)								500		500	
Marijuana (THC)	50					50	50	50	50	50	
Methadone (MTD)								300	300	300	
Methadone Metabolite (EDDP)											
Methamphetamines (MET)						1000	1000	1000	500	1000	
Opiates (MOP) or (OPI300)						300	300	300	300	300	
Oxycodone (OXY)							100	100	100	100	
Phencyclidine (PCP)											25
Tricyclic Antidepressants (TCA)											
** Cotinine / Nicotine (COT)	200										
** Ethyl Glucuronide (ETG)		500									
** Fentanyl (FEN)			25								
** Gabapentin (GAB)				2000							
** Kratom (KRA)					100						
** Synthetic Marijuana (K2)											
** Tramadol (TRA)											

** NOTE: Premium panel not CLIA Waived intended for forensic use only.

A green box  indicates the substance in the far left column is tested.

The number inside the box indicates the cut-off level in ng/ml. A cut-off level is the concentration a med or substance must exceed in the urine sample to trigger a positive result.

Learn more about cut-off levels & detection times in this image gallery.

Tests above with the American Flag  are designed & manufactured in the United States.

All Identify Diagnostics & Identify Diagnostics USA tests are made with USA reagents.

Reagents are the substances on each panel strip which produce the urine test result.

Identify
Diagnostics

ID-CP12-BUP



Identify
Diagnostics USA

ID-US12-BUP




Identify
Health


ID-H12-1



5 PANEL
(single sided)



12 PANEL
(double sided)



1 PANEL
(FEN)
fentanyl



1 PANEL
(ETG)
alcohol



1 PANEL
(COT)
nicotine



10 PANEL
(double sided)



Identify
Diagnostics

Identify
Diagnostics USA

Target Markets

In Office Physician Labs- CLIA
Employment Testing
Parental / Family Testing
Government
Addiction

Distributed By

Medical Distribution Group, Inc.
2361 Whitfield Park Avenue
Sarasota, FL 34243
P: 727-744-2967
F: 727-499-7355

HEP-C



HIV



LC-MS SUPPLY



CORFLEX BRACES



CareFull catch BY CAREFULL PRODUCTS

