

Identify[®] CUP TESTS

* Forensic Use Only Panel or Device.

ID-CP5 ID-CP6 ID-H6-1 ID-H7-1 ID-H9-1 ID-CP10 ID-CP12-BUP ID-CP12-BUP-ADULT ID-US12-BUP MADE IN USA ID-H12-EFKK-1 * ID-H12-1 ID-H12-2 ID-H13-1 ID-US14-ADULT MADE IN USA ID-H14-ADULT-1 ID-H14-ADULT-EFKT-2 * ID-H16-ADULT-EFKT *

SUBSTANCE TESTED	CUP TESTS																
Amphetamines (AMP)	1000	500	500	500	1000	1000	1000	500	500	500	500	500	500	500	500	500	500
Barbiturates (BAR)				300	300	300	300	300		300	300	300	300	300		300	
Benzodiazepines (BZO)	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300
Buprenorphine (BUP)			10	10		10	10	10	10	10	10	10	10	10	10	10	10
Cocaine (COC)	300	300	150	150	150	300	300	300	150	150	150	150	150	150	150	150	150
Ecstasy (MDMA)						500	500	500		500	500	500	500	500	500	500	500
Marijuana (THC)	50	50	50	50		50	50	50	50	50	50	50	50	50	50	50	50
Methadone (MTD)				300	300	300	300	300		300	300	300	300	300	300	300	
Methadone Metabolite (EDDP)													300	300		300	
Methamphetamines (MET)				500	1000	1000	1000	500	500	500	500	500	500	500	500	500	500
Opiates (MOP) or (OPI300)	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300
Oxycodone (OXY)	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100
Phencyclidine (PCP)					25	25	25	25		25		25	25	25		25	
Tricyclic Antidepressants (TCA)											1000	1000	1000	1000			
* Ethyl Glucuronide (ETG)									500							500	500
* Fentanyl (FEN)									25							20	25
* Kratom (KRA)									100								
* Synthetic Marijuana (K2)									25							25	25
* Tramadol (TRA)																100	100

ADULTERANTS

Creatinine (C)									✓						✓	✓	✓	✓
Glutaraldehyde (G)															✓	✓		✓
Acidic or Alkaline pH (P)									✓						✓	✓	✓	✓
Nitrite (N)															✓	✓		✓
Oxidants (O)															✓	✓	✓	✓
Specific Gravity (S)									✓						✓	✓		✓

12 PANEL



Identify[®]
Diagnostics

14 PANEL



Identify[®]
Diagnostics USA

14 PANEL





Identify[®]
Health

Identify[®] DIP TESTS

* Forensic Use Only
Panel or Device.

ID-H1-COT-DIP *
ID-H1-ETG-DIP *
ID-H1-FEN-DIP *
ID-H1-GAB-DIP *
ID-H1-KRA-DIP *
ID-CP5-DIP
ID-CP10-DIP
ID-H10-DIP
ID-CP12-DIP
ID-H12-EFT-DIP *

SUBSTANCE TESTED	DIP TESTS									
Amphetamines (AMP)					1000	1000	500	1000	500	
Barbiturates (BAR)						300	300	300		
Benzodiazepines (BZO)						300	300	300	300	
Buprenorphine (BUP)								10	10	10
Cocaine (COC)						300	300	150	300	150
Ecstasy (MDMA)							500		500	
Marijuana (THC)						50	50	50	50	50
Methadone (MTD)							300	300	300	300
Methamphetamines (MET)						1000	1000	500	1000	500
Opiates (MOP) or (OPI300)						300	300	300	300	300
Oxycodone (OXY)							100	100	100	100
Phencyclidine (PCP)									25	
* Cotinine / Nicotine (COT)	200									
* Ethyl Glucuronide (ETG)		500								500
* Fentanyl (FEN)			25							50
* Gabapentin (GAB)				2000						
* Kratom (KRA)					100					
* Tramadol (TRA)										100



1 PANEL
(FEN)
fentanyl



1 PANEL
(ETG)
alcohol



1 PANEL
(GAB)
gabapentin



12 PANEL
double sided
(ETG) (FEN) (TRA)
alcohol fentanyl tramadol



* Forensic Use Only Panel or Device.

DT-AT12-1

DT-AT14-AD6-1

DT-AT18-FLAT-1

DT-AT20 *

SUBSTANCE TESTED	CUP TESTS			
	DT-AT12-1	DT-AT14-AD6-1	DT-AT18-FLAT-1	DT-AT20 *
Amphetamines (AMP)	500	500	500	500
Barbiturates (BAR)	300	300	300	200
Benzodiazepines (BZO)	300	300	300	300
Buprenorphine (BUP)	10	10	10	10
Cocaine (COC)	150	150	150	150
Ecstasy (MDMA)	500	500	500	500
* Ethyl Glucuronide / Alcohol (ETG)				300
Fentanyl (FYL)	1	1	1	100
* Gabapentin (GAB)				2000
Heroin-Metabolite (6-MAM)			10	100
* Ketamine (KET)				1000
* Kratom (KRA)				300
Marijuana (THC)	50	50	50	50
Methadone (MTD)	300	300	300	300
Methadone Metabolite (EDDP)			300	
Methamphetamines (MET)	500	500	500	500
Morphine (MOP)	300	300	300	300
Oxycodone (OXY)	100	100	100	100
Phencyclidine (PCP)		25	25	25
Propoxyphene (PPX)			300	
* Synthetic Marijuana (K2)				50
Tramadol (TML)			100	200
Tricyclic Antidepressants (TCA)		1000	1000	

ADULTERANTS

Creatinine (C)	✓		✓
Glutaraldehyde (G)	✓		✓
Acidic or Alkaline pH (P)	✓		✓
Nitrite (N)	✓		✓
Oxidants (O)	✓		✓
Specific Gravity (S)	✓		✓

12, 14, 20 PANEL ROUND CUP



18 PANEL FLAT FACE CUP

